

## REAL-TIME MANAGEMENT FROM DOSE TO INVOICE

## Radiopharmacy Management Information System



- Ordering and Billing Solution
- Inventory Management
- Real time Dispensing
- Task List Management

**RMIS** is a management system designed to fit the needs of nuclear pharmacies.

It will streamline inventory receipts, processing prescription orders, shipping, and invoicing. **RMIS** will allow you to efficiently operate your nuclear pharmacy and to comply to your local regulations.


Reduce costly transcription errors and minimize stress for your staff with our RMIS.

A simple solution for nuclear pharmacies. RMIS streamlines your business processes and enhances customer service.

The software provides an electronic ordering and dosage information system to your customers. It combines a stock management, real-time dispensing monitoring, a dose information system and a flexible billing module that integrates with OuickBooks.

The software can eliminate delivery errors with GPS tracking, transcription errors with an online ordering option, and any other errors with a dose calibrator quality control system.

The digital task list and the centralized data enable pharmacy staff to efficiently process orders, track inventory and manage health physics.



Official partner of:







## Features are:

- Direct communication with dose calibrators and real-time dispensing documentation complete with precise filled amounts.
- A calculator to determine the amounts needed to fill orders while providing the ability to monitor activity and adjust it quickly.
- User-friendly interface that includes touch screens enclosed in hoods, high-speed printers and bi-directional communication with dose calibrator.
- Print syringe and Rx labels
- Quick access to an elute generator and mix kits screen
- Stock management to view inventory items, all on one screen, track expired products and other inventory details, print labels and easily transfer products between locations.
- A tool for barcode scanning, kit preparation, and quality control.
- Integration with barcode scanners inside hoods, and at shipping stations to verify correct information input.
- A centralized place to monitor dose calibrator constancy, linearity and accuracy, sealed source inventory, survey meter calibrations, area monitors and task lists.
- Reporting capabilities for kit preparation, dispensing, billing, quality control, inventory and shipping.
- Detailed reports and compliance with 797 Live Dispensed Activity (LDA) regulations, always ensuring precise inventory and RMIS validation.
- An export module to QuickBooks for effortless billing and recordkeeping.
- A flexible billing module allowing to create and email invoices, bill by customer or by cycle, and apply service charges or credits to accounts as needed.
- Backup Redundance by 3 diverse ways of backing up data. (Raid hard drive, Mirroring, and redundant data base server).

Hospitals with NMIS can place their orders online, which are sent directly to the RMIS program at the pharmacy. You'll receive order alerts, and customers receive updates right on their scheduler.

For more details about our products and services, please contact us.



info@elysia-raytest.com www.elysia-raytest.com



Headquarters

Elysia s.a. Rue du Sart-Tilman 375 4031 Angleur - Belgium Tel +32 (0)4 243 43 50

Production Elysia-raytest GmbH Benzstraße 4 75334 Straubenhardt - Germany Tel +49 (0)7082 92 55 0